



The **Medtronic ACT II** is an Automated Coagulation timer. The Act II is a portable device making it easy to run coagulation test in an OR while performing a cardiovascular or vascular surgery, as well as catheterization in laboratories and in Critical Care Units. The with Medtronic ACT II you can run a variety of test including Activated Clotting Time (ACT), High Range Heparinase (HRHTC), Prothrombin Time (PT), a General Purpose Cartridge (GPC) and the Activated Partial Thromboplastin (aPTT). The Coagulation Timer come with a DPU-414 thermal printer.



Specifications

Dimensions

Height: 6.5 in (16.5 cm)
Width: 6.5 in (16.5 cm)
Depth: 9.5 in (24.1 cm)
Weight: 8lbs (3.6 kg)

Other Attributes

Citrated whole blood coagulation controls for verifying instrument and cartridge performance; flashing digital display and audible signal indicate end of test; automated incubation is standard; heparinase cartridge measures ACT of sample with and without

Display

Digital LED

Used (Staining Posts)

Dialysis, surgery, cardiac cath lab, physician lab, clinical lab, ICU

