



The MEDTRONIC BIO-MEDICUS Model 550 Bio-Console Pump Speed Controller offers control and monitoring for the Bio-Pump Centrifugal Blood Pump. Ergonomic design features uncluttered control panel and highly visible graphics. Compatible with most heart-lung consoles. RPM control knob alerts operator that a low-flow condition is approaching.

Features

- Analog and Digital displays for flow and RPM
- Color LCD display for greater ease of use
- Now with TWO pressure measurement capabilities
- New distinctive user choices for alerts/alarms
- Updated flow control knob for enhanced ergonomics
- Remote capability for convenient LCD/flow controller placement
- Unique base unit display offers redundancy for peace of mind
- Updated battery charge status display
- Updated electronic architecture provides a bright future for new functional enhancements



Specifications

General

Medtronic Bio-Medicus Bioconsole

Centrifugal Extracorporeal Blood Pump

Model: 550

Serial No: 7530

120 VAC – 60 Hz – 3A

Controls: Pump Speed, 2 Timers, Alarm Limits, Pressure zero, Volume

Switches: Master Power, Display Test, Motor Transfer, AC/28V/Battery, Battery Test/Standby

Displays: Pump Speed, Flow, Pressure, Battery, 2 Timers

Hospital grade cord

Red protective pump cap