



The **Covidien Nellcor WarmTouch 5300A** is a convective air warming system that keeps patients warm and minimizes postsurgical complications and recovery times. Comfortable, durable blankets and quiet, efficient warming unit give clinicians a convenient, cost-effective way to maintain patient temperatures before, during, and after surgery. It has multiple mounting options including on a cart, bedrail, IV pole, and the floor.

## Features

- Warms up to 38°C in just 30 seconds.
- Boost setting for patients who need fast, short-term warming.
- High airflow of 55 CFM.
- HEPA filter to reduce airborne particulates.
- Comfortable and durable blankets.



## Specifications

### Dimensions

**Height:** 15" (38 cm)  
**Width:** 16" (41 cm)  
**Depth:** 11" (28 cm)  
**Weight:** 15 lbs. (7.3 kg)

### General

**Noise Level:** Less than 55 db  
**Airflow Range:** About 55 CFM  
**Warm-Up Time:** Air temperature nozzle exceeds 38°C in about 30 seconds.  
**Heater Size:** 1,000 watts  
**Hose:** 5 ft. (1.5 m); Collapsible and self-supporting

### Blankets

**Adult**  
Upper Body  
Lower Body  
Full Body

**Pediatric**  
Small  
Large

**Sterile Cardiac**  
Cardiac

