



The **Medtronic Valleylab™ FT10 Energy Platform** not only drives our full portfolio of energy-based devices, it makes them better than ever. Improving the performance of all of our devices with faster sealing and cutting times — for example, LigaSure™ vessel sealing is up to 50% faster — the Valleylab™ FT10 energy platform delivers precise amounts of energy where and when it is needed most.

Features

- 50% faster sealing times
- Lower jaw temperatures
- Auto bipolar has a faster activation time
- Monopolar performance is improved
- Power by smart generators that manage every energy delivery
- Improves all surgical applications



Specifications

Physical Specifications

Weight (kg): 10.1

Width (cm): 36.8

Height (cm): 17.8

Depth (cm): 46.2

Performance

Lower max jaw temperature

Faster sealing times: (1-4 seconds vs. 3-6 seconds with Force-Triad™)

Automatic power settings require minimal setup and minimize need for further handling during surgery

Simple, intuitive controls and information displays

New ES Settings

Soft coag

TissueFect

Reads tissue: 434,000 times per second

Bipolar Cable Compensation

Reads cable length and width for consistent electrosurgical output

