



The **Medtronic Midas Rex Legend EC200** is an electric high-speed drill. The Midas Rex drill is one of the most popular systems used in spinal procedures. Its small and lightweight design makes it suitable for other procedures including neurotology and ENT.

Features

- Small and Lightweight
- Quick Setup
- Easy-to-use Design
- Adjustable speed
- Powered by the Integrated Power Console



Specifications

Dimensions (Console)

Height: 4 In (10.16 cm)

Width: 10 In (25.4 cm)

Depth: 13 In (33.02 cm)

Power

VAC: 100 - 240

Hz: 50/60

Amps: 2

Speeds

EHS Speed: 10,000 - 70,000 RPM (Increments of 1,000)

Stylus Speed: 3,000 - 9,000 RPM (Increments of 500)

Accessories

1 EA300 Motor Cable

1 EM100A EHS Motor

1 Footswitch

1 Footswitch Cable

3 Attachments

