



The **Midmark 417 Podiatry Chair** has everything you'd want in a chair custom-built for your practice: contemporary styling with clean, handsome lines - and feature-packed convenience for podiatrist, patient and staff! Its foot control couldn't be easier to program, or more responsive. A touch is all it takes to position each patient's feet precisely where you want them. Or re-position them during procedures without compromising your sterile field. And it won't slip, slide - or scuff the floor. The optional Hand-held Commander mounts on either side of the foot section to provide 7 additional choices.

## Features

- Adjustable height, tilt, - work standing or sitting.
- Fastest power podiatric chair - 1 button control.
- Foot control (standard) - 3 power positioning, 2 programmed choices.
- Hand-Held Option - 3 power positioning, 4 programmed choices.
- One-Touch automatic return - saves seconds on each procedure.
- Deep, sculptured cushions - built-in pillow.
- Duplex hospital-grade electrical receptacle - left and right.
- Vinyl upholstery - stain-resistant; Velcro-fastened.
- ISO 9001 certified.



# Specifications

**Minimum Height**

**Chair Position:** 22-1/2"  
**Table w/o tilt:** 19-1/2"

**Maximum Height**

**Chair position:** 41"  
**Table w/o tilt:** 38-1/2"  
**Table w/maximum tilt:** 38-1/2"  
**Table w/maximum tilt, maximum lift, foot section fully extended:** 58"

**Length**

**Back reclined:** 70"  
**Foot section fully extended:** 78.69"  
**Foot section:** 14" long

**Width**

**Back:** 22-3/4" tapering to 17-1/2"  
**Seat:** 22-3/4" tapering to 21-1/2"  
**Foot:** 21-1/4" tapering to 17-3/4"  
**Total w/arms:** 28.12"  
**Arm rest width:** 3.12"

**Debris Tray**

8-3/4" W x 12-1/2" L x 1-3/4" D  
**Weight capacity:** 300 lb

**Electrical**

115 VAC, 60 Hz, 12 amps  
Two hospital-grade duplex electrical receptacles, rates 115 VAC, 60 Hz, maximum 6 amps.

**Chair Options**

Stainless-steel debris tray, 9-1/2" W x 11-1/2" L x 2-1/2" D  
Vinyl scuff cover for foot rest

