



The **Midmark 623 Barrier-Free Examination Table** allows patients incapable of getting onto an exam table (whether it be the elderly, or young children) to easily get on, increasing the patient's comfort and satisfaction with the patient experience. To meet the needs of America's disabled patients, as well as the physicians who care for them, the barrier-free exam tables decrease the chance of injury, helping maintain the patient's dignity and reducing their anxiety. This will help increase exam room efficiency, comfort, and safety.

Features

- Low-height
- Heating control for upholstery
- Accessible heating controls
- Premium Comfort System material
- Foot control
- Hand control
- Easy-to-clean removable upholstery



Specifications

Physical Specifications

Top

- 28.38" W x 59.5" L (72.1 cm x 151.1 cm)
- Length with Footrest Extended:** 77.5" (196.9 cm)

Height

- **Minimum:** 18" (45.7 cm)
- **Maximum:** 37" (94.0 cm)

Patient Load Rating 400 lbs (181.4 kg)

Paper Roll Holder Holds 21" x 3.5" paper roll (53.3 cm x 8.9 cm)

Rear Storage Pod 17.8" x 7.0" x 9.8" high (45.2 cm x 17.8 cm x 24.8 cm)

Power

115 VAC
50/60 Hz
12.5 Amps
UL 60601-1 and CAN/CSA-C22.2 No.601.1-M90

