



The **Midmark 641** is a procedure chair that offers four axis power positions. This allows for better ergonomics for the medical professional while still providing comfort for the patient. Additionally, the 641 offers unobstructed access to the patient. With its 19-inch low height, it is easy for patients to access the chair with little to no assistance. The procedure chair was designed to naturally be like a chair while still providing head and neck support. The Midmark procedure chair offers a power height adjustment that can go from 19 inches to 40 inches in height. The 641 also provides a power-adjustable back, allowing the procedure chair to transfer to positions from flat, semi-flat, or a Trendelenburg position.

## Features

- Close-up Access: Narrow patient support surface, slim back section and tuck-away side rails
- Premium Positioning: Put the chair exactly where you need it with powered seat, back, foot and tilt
- Comfort: Exclusive Premium Comfort System enhances patient comfort with a chair-like feel
- Programmability (Optional): four positions that you can customize to your needs
- Electrical Outlets (Optional): duplex, hospital grade (on each side)
- Rotation (Optional): allows 350 degrees of chair rotation and adds 3.5" to the minimum and maximum seat height specifications
- Wireless Controls (Optional): hand and foot
- Patient Weight Capacity: 450 lb (204 kg)
- Height: 19" (48.3 cm) minimum - 40" (101.6 cm) maximum



# Specifications

## Dimensions

Height: minimum 19 in (48.3 cm)  
Height: maximum 40 in (101.6 cm)  
Length: 73 in (185.42 cm)  
Length: 79 in (200.66 cm) (maximum with headrest extended)  
Width: 23.375 in (59.37 cm)

## Procedure Chair

Patient weight Capacity: 450 lbs (204 kg)  
Seat Tilt: 0° to 30°  
Back support angle range: 0° to 85°  
Foot Support Angle Range: 0° to 90°  
Rotation (optional): allows 350° of chair rotation  
Upholstery Width Dimensions:  
Headrest: 9 in (22.9 cm)  
Back support: 21.5 in (54.6 cm)  
Seat: 23.38 in (59.4 cm)  
Foot support: 21.75 in (55.2 cm)

## Electrical

115/230 VAC, 8.5/4.5 A, 50/60 Hz (for standard and rotation models)  
115 VAC, 11.5 A, 50/60 Hz (for standard models w/outlets)  
Optional Electrical Outlet Rating: 115 VAC, 3 A

