



The **Midmark M11 UltraClave** is self-contained steam sterilizer that can be used in medical and dental offices, hospitals, clinics, nursing homes, laboratories, and other facilities to sterilize heat and moisture stable reusable items. The M11 UltraClave is designed to optimize sterilization output, reducing time and effort.

Features

- Dry Cycle Door Stop - Prevents door from opening fully if there is a slight residual pressure in the chamber.
- Monitors chamber temperature during a cycle to prevent an overheat condition.
- Pressure relief valve opens to reduce chamber pressure if the chamber steam pressure exceeds allowable limit.



Specifications

Dimensions

Height (With Printer): 17.9 in. (45.2 cm)
Width: 17.8 in. (45.2 cm)
Depth (With Plug): 23.8 in. (60.5 cm)
Counter Area: 17.8 x 21 in (45.2 x 53.3 cm)
Chamber: 11.0 in diam. x 18 in deep (22.9 cm diam. x 46 cm deep)
Large Tray: 9 in x 15 in x 1 1/8 in (22.9 cm x 38 cm x 2.9 cm)
Small Tray: 6 5/8 in x 15 in x 1 1/8 in (14.3 cm x 38 cm x 2.9 cm)
Weight (Empty): 99 lbs (44.9 kg)
Weight (Full): 111 lbs (50.4 kg)

Water Reservoir Capacity

Approximately 1.4 Gallons (5.3 liters) to full mark. Usable volume is 1 Gallon (3.8 liters)

Electrical Requirements

115 VAC, 50/60 Hz, 15 Amps Single Phase
115 VAC has 1425 watt tubular immersion heater
230 VAC, 50/60 Hz, 10 Amps Single Phase
230 VAC has 1500 watt tubular immersion heater

Chamber Pressure

At 270 °F (132 °C) : 27.1 psi. (182 kPa)
Safety Valve Setting: 40 psi. (275.8 kPa)

Heat Emission

5000 BTU/Hr during operation

Fuses

115 VAC models
F1: 0.250 amp, 250 V, Slo-Blo, 1/4" x 1-1/4"
F2: 15 amp, 250 V, Fast-Acting, 1/4" x 1-1/4"
230 VAC models:
F1: 0.125 amp, 250 V, Slo-Blo, 5mm x 20mm
F2: 8 amp, 250 V, Fast-Acting, 5mm x 20mm

