



The **Midmark Ritter 204** offers long term value with its cost efficiency and comfortable patient care. The Midmark Ritter 204 manual exam tables have 3 pass thru side drawers and 2 front drawers for optimal supply management. The exam table has a reinforced steel seat designed to prevent deflection and breakage. The upholstered top is easily removed without tools via patented knob system.

Features

- 500 lbs weight capacity.
- Manual back with spring assist.
- Steel reinforced seat.
- Drawers manufactured with one piece polystyrene construction.
- Ergonomic design.



Specifications

Height with Upholstery **Soft Touch Top:** 32.13" (81.6 cm)
Seamless Top: 33.13" (84.2 cm)

Patient Weight Capacity **Capacity:** 500 lbs (226.8 kg)

Electrical Receptacle Duplex, hospital grade electrical receptacle on left side

Soft Touch Upholstery Top **Width:** 30.5" (77.5 cm)
Length: 56" (142.5 cm)
Paper Roll Holder: Holds 21 x 4.75" Paper Roll (53.3 x 12 cm)
Length With Footrest: 73.75" (187.3 cm)

Seamless Upholstery Top **Width:** 28.5" (72.4 cm)
Length: 59.5" (151.1 cm)
Paper Roll Holder: Holds 21 x 4.75" Paper Roll (53.3 x 12 cm)
Length With Footrest: 76" (193 cm)

