



The **Midmark Ritter 222 Barrier Free** powered exam table is designed to lower to just 18-inches from the floor to reduce risk of injury from lifting patients. The Midmark Ritter 222 has a manual gas spring assist back, compared to the Midmark Ritter 223 which has a fully powered adjustable back. The Ritter 222 examination table makes your patients feel more stable and relaxed while they're waiting for their consultation as their feet rest on the floor.

Features

- Lowers to 18-inches from the floor.
- Maximum height of 37-inches.
- Precise and smooth patient positioning with minimal effort.
- Built-in pillows.
- Built-in stirrups that lock in four lateral positions.
- Foot extension shelves to expand table length.
- Foot control.



Specifications

Dimensions

Height: 18" to 37" (45.7 to 94 cm)
Width: 30.5" (77.5 cm)
Length (Without foot extension): 59.5" (151.1 cm)
Length (With foot extension): 76" (193 cm)
Weight: 325 lbs (147 kg)
Patient Weight Capacity: 400 lbs (181.4 kg)

Range of Motion

Height: 18" to 37" (45.7 to 94 cm)
Back: 0° to 80° ± 5° (Gas spring assist)
Pelvic Tilt: Max: 7°; Min: 1.5°

Paper Roll Holder

Holds 21" x 3.5" paper rolls (53.3 cm x 8.9 cm)

Electrical

115 VAC, 50/60 Hz, 11 Amps (with heater & receptacle -016)
230 VAC, 50/60 Hz, 2.75 Amps (basic unit -017, excludes receptacle)

Safety Standards

UL 60601-1
CAN/CSA No. 601.1-M90

