



The **Midmark Ritter 223 Barrier Free** is a power examination table. The 223 Barrier Free is used as a chair or table to provide positioning and support of patients during general examinations and procedures conducted by qualified medical professionals in an office environment. This table is designed to lower to 18 inches from the floor making it easier for patients to get on the table and eliminate access strain on their backs as well as the back of the physicians and staff.

Features

- Power height and back.
- 400-pound patient weight capacity.
- 18" to 37" seat height adjustment.
- Standard foot control.
- Weight: 315 lbs.



Specifications

Dimensions

Height: 18 to 37" (45.7 to 94 cm)

Width: 28.5" (72.4 cm)

Length: 76" (193 cm)

Weight: 315 lbs. (142.8 kg)

Patient Weight Capacity: 400 lbs. (181.4 kg)

Movements

Height: 10 to 37" (45.7 to 94 cm)

Back: 0° (Horizontal) to + 80° ±5°

Pelvic Tilt: Max: 7° ±2°; Min: 1.5°

Height Speed: 12 - 26 seconds

Back Speed: 15 - 20 seconds

Power

115 VAC, 50/60 Hz, 12.5 Amps (with heater & receptacle -016)

