



The **Midmark Ritter 244 Barrier Free** is a bariatric treatment table. The Midmark Ritter 244 has an outstanding 850-pound patient weight capacity and features an 18-inch low height. This table is designed to safely and effectively get your patients on and off the table with minimal assistance. With a maximum height of 34-inches, the Ritter 244 lifts your patients to a comfortable working height, helping minimize stress and strain on your back.

Features

- 32" wide seat.
- Seamless upholstery looks great and is easy to clean.
- Barrier-Free accessibility for your patients.
- Fully powered height, back, tilt, and angle adjustments.



Specifications

Dimensions

Height: 18 to 34" (45.7 to 86.4 cm)

Width (Cushion): 32" (81.28 cm)

Length: 80" (203.2 cm)

Patient Weight Capacity: 850 lbs. (385 kg)

Movements

Height: 10 to 34" (45.7 to 86.4 cm)

Back: 0° (Horizontal) to + 65°

Seat Tilt: Max: ±15°

Upholstery Dimensions

Seat and Back Area: 32" W (81.3 cm)

Head and Foot Area: 30" W (76.2 cm)

Overall Length: 80" (203.2 cm)

Power

120 VAC, 2.5 A, 60 Hz

