



The **Mindray Datascope Gas Module SE** is a companion to several patient monitors. The Gas Module SE, commonly called the Datascope GM SE, is capable of automatically identifying and measuring five anesthetic agents (Desflurane, Enflurane, Sevoflurane, Halothane, and Isoflurane, as well as N<sub>2</sub>O, CO<sub>2</sub>, and O<sub>2</sub>). This unit monitors all gases via sidestream sampling.

## Features

- 200 mL/min sample flow rate.
- Uses a non-dispersive infrared analyzer measuring absorption of the gas sample at seven infrared wavelengths.
- 2.5 second sampling delay
- Water trap that removes water and impurities from the patient sample.
- Analyzes patient sample at a fixed rate.



