

# Mizuho OSI 5855 Orthopedic Trauma Top Modular Table System



The **Mizuho OSI 5855 Orthopedic Trauma Top** for the Mizuho OSI 5803 Advanced Control Base has a carbon fiber construction and provides easy patient positioning. The Trauma Top with the Advanced Control base creates a complete orthopedic trauma table package. This Jackson Trauma Table provides unique skin and skeletal traction during key trauma procedures.

## Features

- Carbon composite construction adds strength and reduces weight
- Steer casters aid table transport
- Unique composite top enhances radiography
- Allows upper and lower extremity imaging at any angle
- Traction Arc™ provides constant traction on patient's leg for adduction/abduction
- Eliminates need for leg spar on operative side
- Permits unrestricted access for surgeon and C-Arm
- Hand controller facilitates height, trendelenburg / reverse trendelenburg adjustments



# Mizuho OSI 5855 Orthopedic Trauma Top Modular Table System

## Specifications

### Dimensions (Base)

**Height (Adjustable):** 27 to 48 In (69 to 122 cm)

**Width:** 32 In (81 cm)

**Length:** 64 to 102 In (163 to 259 cm)

**Patient Weight Capacity:** 500 lbs (226 kg)

### Power

**Voltage:** 120V

**Frequency:** 60 Hz

### Movements

**Trendelenburg /Rev Trendelenburg:** 10°

**Lateral Tilt (Roll):** 25°

### Dimensions (Table Top)

**Length:** 84 In (213 cm)

**Width:** 21.5 In (55 cm)

