



The refurbished **Mizuho OSI 5979-1 Cervical Management Base Unit** is used in conjunction with Mayfield or DORO skull clamps. It's mounted on the crossbar of the Trios Base or the MTS Advanced Control Base equipped with the Radiolucent Top or the Spinal Top. This Mizuho Base Unit provides traction for the cervical spine with a longitudinal slide and fine vertical adjustments with the crank.

Features

- Mounts to the crossbar.
- Compatible with the Trios and MTS base.
- Compatible with Mayfield and DORO skull clamps.
- Compatible with the Radiolucent Imaging Top and the Spinal Surgery Top.

