



During procedures requiring surgical slush, **ORS - 1075HS Hush Slush Systems** from Ecolab® deliver velvet-soft slush which improves patient safety and saves valuable nursing time. The proprietary Auto-Stir® technology in Hush Slush Systems automatically stirs the solution while cooling and eliminates concerns with conventional slush that may have sharp edges which can cause trauma to delicate tissues.

## Features

- Slush Basin with Auto-Stir
- Disc-Drap Technology
- Optional Ecolab ORS Fluid Warming System
- Minimize the Risk of Trauma to Delicate Tissues
- Save Valuable Nursing Time
- Easy to Operate and Maintain
- Quiet Operation
- Portable



# Specifications

## Physical Specifications

**Height:** 36"  
**Width:** 18"  
**Depth:** 18"

## Fluid Level

1 to 4 Liters

## Basic Capacity

7.5 Liters (Round)

## Rear Storage Pod

17.8" x 7.0" x 9.8" high (45.2 cm x 17.8 cm x 24.8 cm)

