



The **Oakworks CFPM 400 Imaging Table** is a 4 movement imaging table designed to accommodate mobile c-arm and ceiling suspended imaging systems. It is designed to accommodate low height for patient wheelchair transfers and is mounted on a low profile base. This imaging table is secured via 4 locking casters. The CFPM 400 comes in two distinct models, each featuring a unique top configuration. The first model boasts a top equipped with an integrated headrest, providing added comfort and support. Meanwhile, the second model features a simple rectangular table top without the integrated headrest.

Features

- Carbon Fiber Top - 56 in. Metal free imaging space
- Top Dimensions - 22 in. (56cm) wide x 84 in. (214cm) long
- Motion 1 - Electric Height Range 26in. - 44in. (66-112cm)
- Motion 2 - Electric Trendelenburg tilt $+15^{\circ}/-12^{\circ}$
- Motion 3 - Electric Longitudinal travel 10in. (25cm)
- Motion 4 - Electric Lateral tilt $\pm 15^{\circ}$
- Controls - Multi-function foot control
- Weight Capacity - 500 lb (227kg)



Integrated headrest top



Rectangular top