



The OAKWORKS® CFPM100 Imaging Table is our newly released single tower lift system that provides the imaging market an exceptional value for a 500 lbs patient weight capacity table. This table is ideal for pain management practices, surgical centers, and other multi-disciplinary medical facilities and is compatible with all mobile c-arm imaging systems. The 25" – 43" adjustable electric height range promotes safe and secure wheelchair transfer due to its super low range of 25" while accommodating the higher ergonomic working requirements at 43". Other technology advancements designed in this product include a new soft touch table base material. This durable base material will protect and minimize C-Arm scrapes and damage. Built into the base is a state of the art foot control that can either stay in the base for easy access for table height adjustment or the foot control can be removed from the base and placed on either side of the table to fit your needs. Included in every purchase is the 22" x 84" table top pad along with 3 cranial positioning pads for either prone or supine procedures. Cost effective solution for new institutions, expanding facilities, and start up practices who require a diving board style elevating imaging table.

Features

- Single tower lift system
- Weight capacity: 500lbs
- Compatible with all mobile c-arm imaging systems
- Table Movement: Height Adjustment only
- Superior Image Quality for Fluoroscopy Guided Procedures
- Adjustable Electric Height Range
- Soft touch table base material
- Durable base material Protects and Minimize c-arm scrapes and damage
- Built-in/Removable foot control
- 4 Year Warranty



Specifications

| | |
|-------------------------------|---|
| Dimensions | 214x68x124 mm (WxHxD) |
| Weight | approx. 1.4 kg |
| Battery Type | Rechargeable NiMH battery pack, quick and easy exchange without opening device |
| Battery operating time | Approx. 4 h at 100 ml/h. When prompted, display shows battery status in h and min |

