



The Oakworks CFPM301 is a 3-movement imaging table with a dynamic height range of 25" to 43".

Features

- 3 Powered Motions
- Motions: Height, Trendelenburg tilt, Longitudinal travel
- Patient Weight Capacity: 500 lbs (227 kg)
- Industry's widest adjustable height range of 18" (46 cm) - 25" to 43" (63 cm to 109 cm)
- +15/-11° Trendelenburg tilt



Specifications

Motion 1	Powered 25"-43" (64-109cm) Height Range
Motion 2	Powered +15°/-11° Trendelenburg
Motion 3	Powered 10" (25cm) Longitudinal Travel
Hand & Foot Control	All powered motions
Choice of Tabletop:	
- Integrated Headrest	22" (56cm) W x 84" (214cm) L
- Rectangular	24" (61cm) W x 80" (203cm) L
Table Voltage Options	120V/60Hz or 230V/50Hz
Table Capacity	500 lbs. (227 kg) distributed evenly
4" (10 cm.) Locking Casters	Locking / Unlocking Casters for mobile or permanent positioning
Padding	1" (3cm) Comfort Foam™
Crescent Face Pad	Integrated Headrest Only
Patient Restraint Strap	Must be used during all procedures
Upholstery	Terra Touch™
Table Weight	538 lbs. (244 kg)
Safety Listings	FDA, cETLus, CE marked