



Oakworks CFPMB301 Bariatric Rated Imaging Table with Rectangular Top incorporates a reputation for high quality and reliability. The CFPMB301 features 3 Powered Motions: Height, Trendelenburg Tilt, and Longitudinal Travel. Patient Weight Capacity is 750 lbs. The CFPMB301 has the industry's Lowest Adjustable Height Range!

## Features

- Rectangular Top
- Patient Weight Capacity: 750 lbs
- 3 Powered motions - Bariatric Rated
- Industry's Lowest Adjustable height Range of 26.5" to 44.5"
- 1" Radiolucent Table Top Pad
- 56" of Metal Free Imaging Area
- Patient Safety Strap
- Heavy Duty 5" Casters
- 1 Hand Control: All Motions
- Certifications: FDA listed, cELT approved, CE marked



# Specifications

## Dimensions

**Rectangular Bariatric Top:** 24"W x 80"L

**Patient Safety Strap:** 84"L with 24" of Velcro® at each end

**Table Weight:** 475 lbs.

**Casters:** Heavy Duty 5"

**Patient Weight Capacity:** 750 lbs

## Certifications

FDA listed

cETL approved

CE marked

## Motions

**Motion 1:** Height Range 26.5" - 44.5"

**Motion 2:** Trendelenburg Tilt: +15°/11°

**Motion 3:** Lateral Tilt: ±15°

## Controls

1 each: Foot (all movements) and Hand (all movements)

