



The **Oakworks Echocardiography EA** table offers many ergonomic features and powered movements to allow for easy patient positioning. Electric cardio panels open up to 15-inches wide, which is ideal for left or right-handed scanning. The Echocardiography EA has a height range of 19 to 33-inches which is ideal for echocardiograms and stress echo. This imaging table has a patient lift and load capacity of 550 pounds.

Features

- 550 lb. maximum lift capacity.
- 3" mattress padding.
- 30" wide tabletop.
- Hand-control with a graphical interface.



Specifications

Dimensions

2 Section Top

Width: 30" (76cm)
Length: 72" (183cm)

Backrest Length : 34" (86cm) long

Seat Length : 38" (97cm) long

Weight Capacity: 550 lb (250kg)

Product Weight: 359 lb (163kg)

Movements

Height Range: 19" - 33" (48-83cm)

Trendelenburg/Reverse: Tilt $\pm 18^\circ$

Fowler Backrest: 0-75°

Sliding Cardio Panel: Opens 15" (38cm)

Access Panel Opens: 10" wide (25cm) x 22" long (56cm)

General

Padding: 3" (8cm)

Fabric: PVC free

Casters: 4" (10cm) locking
Electric access back support

Electric sliding cardio panel

Paper roll holder

