



The CLV-U40 is a universal light source for use with all OLYMPUS endoscopes. By actuating the Mode Selector, the light source can be set for any scope. From the EVIS 100, 130, 140, and 200 Series videoscopes to the OES fiberscopes and all rigid endoscopes.

Features

300W Xenon Light Source

- Provides approximately 1.5x greater light intensity for clearer, more vivid images to enhance detailed observation.
- (Only in combination with 200 series scopes.)

Optimum Light Density Adjustment for Video Taping

- Automatically selects the ideal sensitivity required. Manual override offers variable adjustment by the endoscopist.

Automatic Exposure Control for Still Photography

- Automatically sets light source at the optimum exposure for still or Polaroid photography. Manual override gives the endoscopist complete control of this function.

Built-In Emergency Lamp

- The CLV-U40 system includes a built-in 150W Halogen emergency lamp which is automatically actuated in the event of Xenon lamp failure.

Malfunction Indicator

- Any malfunction in the emergency system triggers the panel light malfunction indicator.

