



This product is compatible with standard and big chip endoscopes. It allows the user to enhance, freeze and adjust colour settings to provide clear and accurate images from inside the body. Requires clv-u20/u40.

Features

- The processor is fully compatible with the Olympus EVIS universal light source
- Compatible with EVIS Video scope 200 series
- Compatible with EVIS CLV-U20 light source
- Power Requirements: 120V AC, 50/60Hz
- Dimensions: (W) 450 mm x (H) 90 mm x (D) 485 mm
- Weight: 33.10 lbs.

