



The Olympus CV-240 Video Processor is compatible with standard and big chip endoscopes. The CV-240 allows the user to enhance, freeze, and adjust color settings to provide clear and accurate images from inside the body.

## Features

- Compatible with the full range of endoscopes from EVIS 100, 130 and 140 series endoscopes
- Wide screen image functionality with used in conjunction with Q series endoscopes such as GIF-Q140 and CF-Q140L
- Pre-set scope switch functions. Endoscopist has a choice from up to 20 pre-set endoscope switch functions
- Pre-program patient data. Up to 40 files of patient data can be pre-entered into memory to expedite procedure turn around

