



The ESG-300 is a high-end HF generator dedicated to advanced gastrointestinal and respiratory procedures such as POEM or ESD. The ESG-300 can be combined with the APU-300 for argon plasma applications.

## Features

- At the Somlpe touch of a button - A simple, intuitive touch-screen interface means users
- can easily set up the system to meet their personal needs in just a few steps.
- Offers Your Guidance - The smart assistance messaging system provides simple and easily understandable feedback that empowers users to confidently operate the ESG-300.
- With No strings Attached - The operational freedom provided by the wireless foot switch allows safer placement in the theatre. Simple, hands-free toggling of procedure settings –
- without the need to directly interact with the generator – means a more efficient user workflow.
- Talks your language - A choice of 29 display languages means users can understand and operate the system with ease, reducing the likelihood of any confusion.
  - Touch screen
  - Wireless foot switch
  - Personalized settings
  - Argon PLASMA ready (in combination with APU-300)



# Specifications

## Dimensions

Size:(W x D x H)295 x 400 x 115 mm

Weight: 6.61 kg

Volume: 13,570 cm<sup>3</sup>

## Out Put

Bipolar Output: International Standard - 2 - pin plug (22 min)

Bipolar Output: Erbe Standard - 1 - pin coaxial plug

Monopolar Output: International 3 - pin plug(4 mm diameter)

Monopolar Output: Bovie Standard- 1 - pin plug (8mm Diameter)

Monopolar Output: Banana plug - 1 pin plug (4mm diameter)

## Max. PowerOutput

120 Watts

## Argon PLASMA

ForcedArgon

PulsedArgon Slow and Fast

SmartArgon

