



With 1080 effective scanning lines of picture information, the high definition television camera head delivers picture quality that is more than twice as good as conventional TV. Images are smooth and clear, with reproduction so accurate that details and colors remain sharp and life-like, opening the way to a new era in endoscopy.

Features

- Autoclavability - The camera head is tightly sealed preventing the lens from fogging even after autoclave steam sterilization.
- One Touch Zoom
- Narrow Band Imaging - Enhances the visibility of vessels and other tissues on the mucosal surface.
- Remote Control of Frequently Used Functions - Three switches are built-in, allowing programming of frequently used functions such as electronic zoom and image capture.

Model Number: OTV-S7PROH-HD-12E



Specifications

Size

Camera Head:

- ø41 x 105 mm (from mount surface)

Cable:

- ø6.8 mm x 4 m

Weight

Camera Head:

- 215 g (excluding cable)

Observation

Pickup System:

- Interline type CCD solid-state image pickup

**Classification
(Electro-Magnetic Equipment)**

Degree of Protection Against Electric Shock:

- Type CF

Degree of Protection Against Explosion:

- The product should be kept away from flammable gases

**Cleaning / Disinfection
Sterilization**

Cleaning / Disinfection:

- Immersible in disinfectant solution

Sterilization:

- Ethylene oxide gas sterilization and steam (autoclavable) sterilization possible

**Number of Remote Control
Switches**

3

Model Number

OTV-S7PROH-HD-12E

