



HDTV1080i combines vivid visualization and comfortable control. With 1080 scanning lines of picture information at your fingertips, HDTV1080i dramatically expands the clinical possibilities in urology and gynecology. The images are smooth and clear, with reproduction so accurate that details and colors remain sharp and true to life. At 60g, this is the most lightweight L-shaped HD camera head available, which enables you to use it for long periods of reduced stress and fatigue during HD minimal invasive surgery, e.g. cystoscopy, hysteroscopy, and resection.

Features

- **High resolution** - The advanced video chip delivers unsurpassed detail and color reproduction in HDTV1080 quality during HD minimal invasive surgery.
- **Lightweight body** - HDTV1080i is the lightest URO HD camera head available, which reduces operator stress and fatigue, even during longer procedures.
- **Narrow Band Imaging** - NBI improves visibility of vessels and the muscle surface allowing detection of cancer and lesions.
- **Standalone remote control switch** - The distantly-positioned remote control switch enables easy control of any three functions and reduces the risk of camera shake.

Model Number: OTV-S7PROH-HD-L08E



Specifications

Size

Camera Head:

- Ø21 x 83 mm

Cable:

- Ø3.3 mm x 4 m

Weight

Camera Head:

- 60 g (head part) 20 g (remote control part)

Observation

Pickup System:

- Interline type CCD solid-state image pickup

Classification (Electro-Magnetic Equipment)

Degree of Protection Against Electric Shock:

- Type BF

Cleaning / Disinfection Sterilization

Sterilization:

- Ethylene oxide gas sterilization possible

Model Number

OTV-S7PROH-HD-L08E

