



The Olympus OEP-3 color A5 printer is designed for endoscopy image printing and features a 300 DPI printing resolution. The Olympus OEP-3 also features over 16 million colors and has a compact footprint with low noise. The Olympus OEP-3 color A5 printer is the ideal solution for demanding print applications.

Features

- Quest International offers sales, repair, service, or exchange for this Olympus printer
- For immediate repair service or consultation on exchange or new options, please complete the form to the right
- High resolution 300 dots per inch (DPI)
- Prints 1, 2, 4, 8 or 16 images per page
- Color Resolution: Over 16 million
- Print Resolution: 300 dpi



Specifications

Dimensions	(h x W x D) cm (in) 49 x 37 x 12.5 (19.3 x 14.6 x 4.9)
Weight	kg (lb) 14 (30.9)
Processing Technique	Dry
Power Needed	100-240 VAC, 50/60 Hz
Dicom compliant (maximum Display)	Yes
Yield, Films/Hr	180
Ram, MB (Storage)	8 (8-frame memory)
Interface Options	RS232C (D-sub 25 pin)
Input Ports	1 RGB, 1 com
Pixels (maximum)	1944 x 1022
Spatial Resolution, Pixels/m	295 dpi
Multiple Originals (Films)	9
Contrast Resolution, gray-scale	16.7 million colors
Internal Modulation, Bits	8
Keyboard Functions	Memory page, source/ memory, capture, print, menu, print quantity, clear, display
Format (Power Needed)	1, 2, 4, 6, 8, 16, split print selection
Film Used	Olympus, Sony
Film size, cm (in) (Film)	14.9 x 21.3 (5.9 x 8.4)
Calibration Technique	Automatic
Modalities (Film)	Endoscopy
FDA Clearance (Interference Compensation)	Yes
ME Mark (MDD) (Interference Compensation)	Yes