



The **Pedigo Midmark 540 stretcher** can hold patients up to 350 pounds. The Midmark 540 features hands-free hydraulic height adjustment, 8-inch casters that can lock from any corner, and four-way positioning for IV poles. The stretcher is versatile and is used for ED, OR, PACU, or any medical department that needs a quality stretcher.

Features

- An optional IV caddy gives you hands-free IV pole transport
- An optional IV folding pole that tucks out of the way when not needed and stays connected to the stretcher when in use
- Precision caster wheels make stretches easy to maneuver, even on carpet. Wheels lock from any corner. 5th wheel steering is also available as an additional option
- Pneumatic assist Trendelenburg and Reverse Trendelenburg are standard
- Optional Fowler X-ray capability available separately
- Hands-free hydraulic height adjustment is smooth and easy
- Tuck away Side Rails let you position the 540 flush with hospital beds for zero transfer gap and easier patient transfer
- Mounted at the foot end of the stretcher, the unique Hand Control Trendelenburg with pneumatic assist easily moves from horizontal to positive or reverse positions
- Sturdy solid steel column with single hydraulic cylinder design uses sixteen precision ball-bearings to smoothly raise and lower the patient



Specifications

Dimensions

Length: 83 in (210.8 cm)
Width: 29.5 in (74.9 cm)
Height (Adjustable): 24 in to 33 in (60.9 cm to 83.8 cm)
Mattress: 77 in x 27 in (195.6 cm x 68.5 cm)
Weight
Normal (depending on options): 300 lbs (136.0 kg)
with shipping Cartoon: 404 lbs (183.2 kg)

Casters

8 in Diameter (20.3 cm)
Central locking of the brakes and steering may be accomplished from all four corners of the stretcher.

Fowler

.0° to 90°, pneumatic assisted

Trendelenburg

12°, pneumatic assisted

Reverse Trendeleburg

7°, pneumatic assisted

