



The fully-functional 550/555's versatility and exceptional quality make it the right stretcher for ED, OR, PACU or virtually any department in your facility.

Features

- 500 lbs. weight capacity
- Stainless steel glide and tuck side rails
- Foot End hydraulics
- Single column hydraulic cylinder design
- Pneumatic assist positive and negative Trendelenburg
- 8" casters with central lock and/or steering with controls from all four corners
- Four corner IV pole receptacles
- Full PVC perimeter bumper (grey)
- Epoxy powder coat Finish
- Pneumatic assist Fowler backrest
- Restraint Strap locations:
 - Fowler.... Four 5" x 7/8" docking ports
 - Hip..... Four 35/8" x 7/8" docking ports
 - Leg..... Two 8" x 7/8" docking ports
- True Directional Steering.



Specifications

Technical Specs

Versatile, Full-Functioning Stretcher that are Utilized in the OR, ER, PACU, Transport and Urgent Care Clinics

Trendelenburg (12 Degrees) and Reverse Trendelenburg (-7 Degrees) Positioning, Pneumatic Assist

Hydraulic Height Adjustment (Surface of Stretcher 24" to 33") Pedal is on the foot end of the stretcher

Hydraulic Height Adjustment

Foot pedal is on the foot end of the stretcher

Fold Away Side Rails

Surround Bumper

Central Brake / Steering

4 Wheel Locking

O2 Tank Holder

Dimensions: 33.5" W X 83" L

Patient Capacity is 500 lbs

New Stretcher Pad