



The Magna-Clave features a series of safety measures, like a Pelton & Crane-patented door lock that seals the entirety of the sterilizer's rim, as well as clarity features like a series of indicator lights ("Power On," "Heat On," "Sterile," and "Open Door"). The Magna-Clave sterilizer is a manually-operated sterilizer with a 15-inch internal-diameter by 30-inch-deep, round, stainless-steel chamber. Due to its large capacity, ease of use, and quick cycle speed (completes a flash cycle in 18 minutes from a warm start), it should come as no surprise that Pelton & Crane's Magna-Clave is an industry-leading sterilizing unit.

Features

- Floor autoclave with steam generator
- Compact size
- Large chamber capacity
- Safety lock system
- Simple and easy locking mechanism
- Easy cleaning



Specifications

Dimensions

Width 28" or 71.1 cm
Depth 34.5" or 87.6 cm
Height With Stand 53" or 134.6 cm
Height Without Stand 23" or 58.4 cm

Chamber Dimensions

15" dia x 26" deep (38.7 cm x 66 cm)

Minimum Space Required

With Stand 34" to 40" wide x 56" deep x 65" high
Without Stand 34" to 40" wide x 56" deep x 35" high

Flash Cycle

18 minutes from a warm start and 35 minutes from cold start.

Indicator

Power On, Heat On, Sterile, Open Door and Water Level

Trays

Small: 10.5" x 26.5" x 2".
Large: 15.25" x 26.5" x 6"

Power

Voltage: 208 - 240 VAC.
Power Consumption: 5000W
Current Rating: 30-amp electrical circuit

