

Penlon Sigma Delta Sevoflurane Anesthesia Vaporizer



The **Penlon Sigma Delta Sevoflurane** is designed for incorporation in the fresh gas supply system of continuous flow anesthesia machines, directly connected between the flowmeter unit and the common gas outlet of the machine. Its purpose is the provision of accurate concentrations of Sevoflurane in the fresh gas supply.

Features

- Weight: 4.8 kg / 10.6 lbs
- Dimensions (HxWxD): 24.2cm x 12cm x 19cm
- Liquid Capacity: 250 mL
- Flow Range: 0.2 to 15 L/min
- Concentration Range: 0.2 to 8 vol.%
- Filling System: Keyed, Screw Top, Quik Fil Filler

