



The **Penlon Sigma Elite Sevoflurane** is an anesthesia vaporizer that is designed for incorporation in the fresh gas supply system of continuous flow anesthesia machines. The Sigma Elite is connected between the flowmeter unit and the common gas outlet of the machine. Its purpose is the provision of accurate concentrations of anesthetic Sevoflurane in the fresh gas supply.

Features

- Weight: <7.3 kg / <16.1 lbs
- Dimensions (HxWxD): 21.7cm x 12cm x 19.4cm
- Liquid Capacity: 300 mL
- Flow Range: 0.2 to 15 L/min
- Concentration Range: 0.2 to 8 vol.%
- Filling System: Keyed or Screw Top Filler

