



Features

- Processor with Xenon light source
- Employs a high speed rotating color filter system producing the most natural color image.
- On-screen data
- Freeze mode
- Infrared capability
- Built-in auxiliary lamp



Specifications

Dimensions

Height: 220 mm
Width: 495 mm
Depth: 540 mm
Weight: 32 kg

Video Output

RGBS, Y/C, Composite (NTSC)

Power Supply

Voltage: NTSC 110V
Frequency: 50-60 Hz
Input Current: 4.8A

Light Source

Energy: 300W
Lamp Type: Xenon

