

## Philips 12-Lead Patient Interface Module (PIM) for TC30/TC50/TC70 Cardiographs



The **Philips 12-lead Patient Interface Module (PIM)** may be configured to support up to 12 optional leads for adult and pediatric application. The same Patient Interface Module (PIM) is used on the TC30, TC50, and TC70 model cardiographs. The PIM is a hand-held device that connects to the patient data cable.

### Features

- Supports 12 optional leads.
- ECG Button that is used to take ECGs from the bedside.
- OEM Part numbers: 453564078291, 453564094261, 453564150741, 453564150761, 453564196651, 453564659321

