

Philips 16-Lead Patient Interface Module (PIM) for TC50/TC70 Cardiographs



The **Philips 16-lead Patient Interface Module (PIM)** may be configured to support up to 16 optional leads for adult and pediatric application. The same Patient Interface Module (PIM) is used on both the TC70 and TC50 model cardiographs. The PIM is a hand-held device that connects to the patient data cable.

Features

- Supports 16 optional leads.
- ECG Button that is used to take ECGs from the bedside.
- OEM Part numbers: 453564196661, 453564150791, 453564150751, 453564150771

