




The **Philips Expression IP5** allows you to monitor patients remotely from the MR Control Room. The Expression IP5 connects to the Invivo Expression and the Expression MR200 patient monitoring systems. The IP5 is used to transmit and remotely integrate your MR patient monitor data from your MRI unit to your hospital's EMR and anesthesia record-keeping systems. Easily access view, and track vital signs without the need for wireless router networks.

## Features

- 19-inch Widescreen display.
- Output 2 sets of data in HL7 and serial data simultaneously.
- Customize and store multiple setups for different procedures and patient types.



# Specifications

## Dimensions

**Height:** 13.04 in (33.12 cm)  
**Width:** 18.75 in (44.04 cm)  
**Depth:** 6.45 in (16.38 cm)  
**Weight:** 18.3 lbs (8.3 kg)

## Display

**Type:** Color WSXGA liquid crystal display  
**Resolution:** 1440 x 900 pixels  
**Touch Screen:** 5 Wire  
**Screen Size (Diagonal):** 19" (48.3 cm)  
**Aspect Ratio:** 16:9 width  
**Refresh Rate:** ≥ 60 Hz

## Input/output Ports (Supported)

**Ethernet:** 1  
**Serial:** 1  
**USB 2.0:** 6

## Peripherals (Supported)

Wireless barcode scanner  
Wireless keyboard and mouse  
USB strip-chart printer  
USB drive (<4GB)

## Power

**Operating Voltage:** 100 to 240 VAC  
**Frequency:** 50/60 Hz  
**Current:** 0.57 A @ 115VAC  
**Power Consumption, Maximum:** 65 Watts

