

Features

-

Specifications

Dimensions

Size: 90 x 370 x 467 mm 3.54 x 14.5 x 18.4 in

Weight: 6.3 kg 13.9 lb

Gas analyzer: M1026BAGM

Power

Power Consumption: peak: 35W, typical: 25W

Power Range: 100 - 240 VAC

Environmental Specifications

Temperature Range

Operating 10 to 40°C (50 to 104°F)

Non-operating -20 to 65°C (-4 to 149°F)

Humidity Range

Operating 5 to 90%

Relative Humidity (RH) max. @ 40°C (104°F) (non-condensing)

Non-operating 5 to 95% Relative Humidity (RH) max. @ 65°C (150°F)

Altitude Range

Operating -305 m to 2900 m (-1000 to 9515 ft)

Non-operating -305 m to 5000 m (-1000 to 16404 ft)

Atmospheric Pressure Range

Operating 70 kPa to 106 kPa

Non-operating 50kPa to 106 kPa

Warmup Time After switching on: 1 - 2 minutes to measure, 6 minutes for full accuracy