



The **Philips IntelliVue MMS Module M3001A** is a module that can be attached to a Patient Monitor. The MMS Module M3001A can provide readings for ECG, SpO₂, NIBP, and temperature. The Philips Module can be connected to the MP or the MX Series of patient monitors via the back of the monitor or a module rack. Making it easy to switch the module from a transport monitor to a bedside monitor without having to adjust the settings. Patient data can be saved on the module and transferred to another monitor. It can store up to 8 hours of patient data on the module.

Specifications

Dimensions Height: 3.75 in (9.525 cm)
 Leigth: 9 in (22.86 cm)
 Width: 2.375 in (6.03 cm)

