



The **Philips IntelliVue MP5CS** is the newest generation in monitoring. The MP5SC provides the medical staff with data that can identify patients at risk. The IntelliVue comes with Guardian Early Warning Score (EWS) that can calculate early warnings and provide a pictorial clue to help respond faster in certain situations. The Monitor goes well beyond the capabilities of a traditional spot-check monitor. what makes this monitor so innovative is the clinical decision support tools that can help measure, evaluate, and document all of the patients' vital signs.

Features

- Essential measurements to support informed decisions
- Large 8.4" touchscreen shows key data at a glance
- IntelliVue user interface simplifies monitoring
- IntelliVue Guardian Early Warning Scoring to aid rapid response
- Flexible operational modes for use on the general ward
- Flexible MX40 connection with short range radio
- Standard supplies – same as rest of IntelliVue monitors
- Clinical Decision Support – Standard IntelliVue Guardian EWS



Specifications

Measuring to support informed decisions

Philips is committed to providing excellent standard measurements such as NBP, Microstream CO₂ (including integrated Pulmonary index), choice of spot check temperature technology with Suretemp, Tympanic Temperature or Temporal Temperature and a choice of pulse oximetry with Philips FAST SpO₂, Nellcor™ OxiMax™ SpO₂ or Masimo® Rainbow® SET® SpO₂ (including certain Masimo rainbow parameters).

Large Touch Screen Shows key data at a glance

This monitor has a large 8.4" touch screen with flexible parameter configurations to enhance patient care.

IntelliVue user interface

The IntelliVue user interface is optimized for ease of use on the general floor. It supports efficient workflow, helps caregivers in a variety of settings to simply measure, evaluate, and document their patient's vital signs and EWS scores.

IntelliVue guardian early warning scoring to rapid aid response

The IntelliVue MP5SC is fully configurable with IntelliVue Guardian Early Warning Scoring (EWS) settings, parameters, alarms, and screen layouts. It provides caregivers with direct access to an automated scoring system based on their own hospital's policy and guidelines. Whatever the level of patient acuity and frequency of vital sign collection, IntelliVue Guardian EWS can help caregivers identify patient deterioration early with fast and easy Early Warning Scoring calculations.

Flexible Operational modes for use on the general ward

The monitor can be used in spot-check or continuous monitoring mode. In continuous mode, the IntelliVue MP5SC can be used for long-term bedside monitoring. Its optional EWS Escalation package allows the MP5SC to convert to a full ECG monitor via a patient-worn device when enhanced monitoring is required, such as during a code event.

Flexible MX40 connection with short range radio

Connects to the IntelliVue MX40 via Short Range Radio for multi-patient vital signs and occasional viewing of continuous monitoring.

Standard supplies - same as rest of IntelliVue monitors

The MP5SC uses the same accessories as the rest of the IntelliVue family, standardizing training, ordering, and use throughout your facility.

Clinical decision support-standard IntelliVue Guardian EWS

Standard IntelliVue Guardian EWS Clinical Decision Support with optional IntelliVue Guardian EWS Escalation Package for patient-worn, device-based wireless ECG monitoring.

