



The **Philips X7-2t TEE Ultrasound Transducer** is compatible with Philips iE33 and CX50 ultrasound systems. The X7-2t probe can be utilized for TEE applications for patients less than 30 kg/66 lb. It features harmonic imaging, as well as dynamic elevational focusing. The array is in excellent condition and can be purchased from Soma Technology, Inc.

## Features

- 2D and 3D Matrix Array
- 2,500 elements
- 7 to 2 MHz
- Harmonic imaging
- Dynamic elevational focusing
- Electrocautery suppression
- Electronic rotatable array from 0 to 180 degrees
- TEE applications for patients >30 kg / 66 lb

