



The Reliance 444 Single-Chamber Washer/Disinfector, a mechanical washer/disinfector, is equipped with an Eagle® 3000 Stage 3™ control system. The washer/disinfector is designed with seven factory set cycles: Instruments, Gentle, Utensils, Glassware, Plastic Goods, Anesthesia/Respiratory Goods, and Decontamination. Three additional cycles are available for customized programming to meet specific operating requirements. Each preprogrammed cycle is equipped with Pre-Wash, Wash, Rinse, and Thermal Rinse (including Drain Discharge Cool Down) phases. The washer/disinfector is built to seismic design and includes a vented system and interior chamber light. The washer/disinfector offers both manifold and non-manifold washing and a power vertical sliding door, in a single- or double-door configuration.

Features

- Available in steam supplied or electric heat configurations as well as vented or no-vented.
- Single door or double vertical door pass-through configurations
- Automatic detergent, enzyme or lubricant dispensing.

The washer is designed with six (6) adjustable cycles each include Pre-wash, Wash, rinse and thermal rinse:

- Instruments
- Gentle
- Utensils
- Glassware
- Plastic goods
- Rubber goods
- Four (4) additional cycles are available for customized programming.



Specifications

Exterior Dimensions

Width: 42 in (1067mm)
Height: 80in (2032 mm)
Length: 32.625 in (829 mm)

Chamber Load Capacity

Width: 24 in (610 mm)
Height: 26in (660 mm)
Length: 24 in (610 mm)

Load Height

31 in (787 mm) above finished floor

