



The **SciCan STATIM 2000** autoclave is designed to provide a complete sterilization solution for unwrapped and wrapped instruments. The STATIM 2000 provides rapid sterilization balanced with rapid drying. The heat from the sterilization process is used to dry the instruments rapidly using a convection method. This convection drying helps dry a properly loaded STATIM cassette.

Features

- Sterilizes handpieces and instruments for immediate use in as little as 6 minutes.
- Thin stainless steel walls of the STATIM's fully-removable cassette, allow for rapid heating and cooling of the chamber.
- Uses convection heat to dry instruments by utilizing the remaining heat in the system after a sterilization process.



Specifications

Machine Dimensions
Height: 6" (15.3 cm)
Width: 16.25" (41.3 cm)
Depth: 17.75" (45.1 cm)
Weight: 46 lbs (20.9 kg)

Cassette Size (Internal)
Height: 1.5" (3.81 cm)
Width: 7" (17.8 cm)
Depth: 11" (27.9 cm)

Power Consumption
Watts: 1,300 watts

Cycle Times
Unwrapped: 6 Minutes*
Wrapped: 12-15 Minutes*
Rubber & Plastics: 20 Minutes*

*Cycle times vary depending on load size

Reservoir Capacity
Reservoir Capacity: 1.05 Gallons (4 Liters)

