



The **Skytron 1602 Essentia** provides 180° top rotation and achieves comfortable surgical or imaging access – all at an impressive value. An ideal solution for versatile operating rooms and surgical centers, this table offers one-touch positioning options for a diverse range of specialties.

## Features

- 180° top rotation
- 500 lb. lift and articulation capacity
- Battery power
- Positioning options include 90° back-up position
- 27-3/16" – 40-15/16" adjustable table height
- 23° lateral tilt
- 28° Trendelenburg and reverse Trendelenburg



## Specifications

**Size**

Length w/o head rest: 68.5" 1,740mm  
 Length with head rest: 78.75" 2,000mm  
 Width (back / seat plate): 19.625" 500mm  
 Width with side rails: 21.875"

**Weight:**

Max. permissible load: 500 lb (227kg) Lift  
 500 lb (227kg) Articulation

**Adjustments via pendant control:**

Height: 27 3/16 to 40 15/16" 690-1,040mm  
 Lateral tilt (left/right): 23°  
 Trendelenburg: 28°  
 Reverse Trendelenburg: 28°  
 Back up/down: 90°/30°  
 Leg down: 100°  
 Return to level (anatomically correct return to horizontal position from articulated positions)  
 Flex/Reflex: 30°/90°  
 Base braking (lock/unlock)

**Adjustments via manual control:**

Top Rotation: 180°  
 Head Section: 60° up/90° down

