



The 3003 Impulse puts automatic control of the imaging surface right at your fingertips. Featuring joystick control, the table offers 65" of unobstructed imaging. The two-speed, multidirectional top slide takes the work out of manual maneuvers and automatically manipulates the imaging surface.

## Features

- 20" longitudinal and 8" lateral top slide
- 500 lb lift and articulation weight capacity
- Two-speed, multidirectional joystick control
- 28"- 42" adjustable table height
- 83" carbon fiber top
- 65" unobstructed imaging window
- 20° lateral tilt
- 20° Trendelenburg and reverse
- Trendelenburg



## Specifications

### Dimensions

**Length:** 83" (2,100 mm)  
**Width:** 19 5/8" (500 mm)  
**Width (With Side Rails):** 22" (560 mm)

### Adjustments Via Hand Control

**Height:** 28"-42" 710-1065 mm  
**Lateral tilt (left/right):** 20°  
**Trendelenburg:** 20°  
**Reverse Trendelenburg:** 20°

### Joystick Control

**For operation, top button on joystick control must be depressed.**  
**Orientation Button:** Reverse Normal  
**Center Button:** Moves table top to center position  
**Speed:** High Low

### Max Permissible Load

500 lb (227kg) Lift  
500 lb (227kg) Articulation

