




Developed for maximum C-Arm access in procedures where only lateral tilt and Trendelenburg positioning are required, the Elite 3100A surgical table provides 23-1/2 inches of longitudinal top slide. This table is the ideal choice for vascular surgery, hip revisions, ureteroscopies, laser angioplasties, pacemaker insertions, pain management and other C-Arm assisted procedures that can be accomplished quickly and efficiently.

## Features

- 500 lb. lift & articulate capacity.
- 23" top slide (toward foot end).
- X-ray cassette tunnel.
- Precision locator system.
- Removable leg section.
- Line power.
- Specialty articulation functions.
- Durable one-touch pendant control.





## One Touch Pendant Control

- Recessed buttons reduce the possibility of accidental activation.
- Durable, water resistant steel housing.
- Pendant control can be used from any point around the perimeter of the table.
- Brakes automatically lock when any button is pushed for patient and staff safety.
- Can be used in conjunction with optional foot control.

## Specifications

### Dimensions

**Length:** 69" (1750 mm)

**Length (with Headrest):** 82" (2080 mm)

**Width: 19.75" (500 mm)**

**Width (with Side Rails):** 21.75" (550 mm)

**Weight of Operating Table:** 774 lbs (351 kg)

**Maximum Permissible Load (Lift & Articulation):** 500 lbs (227 kg)

### Hand Control Adjustments

**Height:** 29" - 40" (737 - 1010 mm)

**Lateral Tilt (Left / Right):** 20°

**Trendelenburg / Reverse Trend:** 20° / 20°

**Power Kidney Bridge:** Yes

**Back Up / Down:** No

**Leg Up / Down:** No

**Top Slide / Head:** 23.5" (596.9 mm)

**Top Slide / Foot:** 23.5" (596.9 mm)

**Flex/Reflex:** No

**Base Braking (Lock / Unlock):** Yes

**Patient Orientation (Normal / Reverse):** Yes

**Return to Level:** Yes

