



Skytron 3501B EZ Slide is a universally versatile surgical table providing excellent imaging flexibility, clinical versatility and cost effectiveness. EZ Slide provides an impressive range of full body imaging and articulating table functions as required for C-arm imaging orientations, without the use of multiple support accessories or having to reposition the patient. One touch of a button does it all!

Features

- 700 lb. lift, 600 lb. articulate capacity
- 21" top slide
- Removable leg section
- Perfect for Vascular, Neuro, Ortho and Urology
- Return-to-level
- Full articulation capacity
- Full body imaging without the need for special accessories
- Durable back-lit pendant control
- Battery power



Specifications

Dimensions

Length w/o head rest: 69 in. (1750mm)
Length with head rest: 82 in. (2080mm)
Width (back / seat plate): 19.75 in. (500mm)
Width with side rails (back / seat plate): 21.75 in. (550mm)
Weight: 827 lbs. (376kg)
Permissible Load: 700 lbs. (320kg) Lift / 600 lbs. (275kg) Articulation

Adjustments Via Hand Control

Height: 26-45 in. 660-1140mm
Lateral tilt (left/right): 20°
Trendelenburg: 25°
Reverse Trendelenburg: 20°
Power Kidney Bridge: 6 in. Lift 150mm Lift
Back up/down: 90°/40°
Leg up/down: 2°/98°
Return to level (anatomically correct return to horizontal position from articulated positions)
Top slide/head: 7.5 in. 195mm
Top slide/foot: 13.5 in. 335mm
Flex/Reflex: 40°/90°



7.5 in. (195mm) Slide to Head



13.5 in. (335mm) Slide to Foot

