



The **Skytron 3600B** UltraSlide provides 23 inch Top Slide, full table features, and optimal patient positioning for C-arm access. It has a powerful 1,000 lb. lift and 800 lb. articulation. Removable back and leg sections and Auto Beach Chair position provide flexible patient positioning options for every procedure. The UltraSlide 3600B also comes with a durable back-lit, one-touch pendant control that is easy to use and provides high visibility when working in low light environments.

Features

- 1,000 lb. lift, 800 lb. articulation capacity.
- Removable back and leg sections.
- Durable, back-lit hand pendant.
- Electrohydraulic positioning.
- Four point braking system with push button operation including manual override.



Specifications

Dimensions

Height: 24-42.5 in (61-108 cm)
Width (Without siderails): 19.5 in (50 cm)
Width (With siderails): 22 in (56 cm)
Length (Without head rest): 69 in (175 cm)
Length (With head rest): 82 in (208 cm)
Weight: 800 lbs (365 kg)

Maximum Loads

Lift: 1,000 lb.
Articulation: 800 lb.

Movements

Height: 24-42.5 in (61-108 cm)
Lateral tilt (left / right): 30°
Trendelenburg: 30°
Reverse Trendelenburg: 30°
Power Kidney Bridge: 5.75 in lift (14.5 cm)
Back (up / down): 90°/40°
Leg down: 15°/100°
Return to level (anatomically correct return to horizontal position from articulated positions)
Top (slide / head): 9.5 inch 24 cm
Top (slide / foot): 13.5 inch 34 cm
Flex / Reflex: 40°/90°
Base braking (lock / unlock)
Beach Chair

