



Pairing the industry's longest top slide with substantial weight capacities, imaging access is only part of the 3602 UltraSlide advantage. Combine the table's top slide capabilities with an 800 pound articulation capacity, and the result is the best of both worlds: imaging and power. The table's impressive 23" of top slide accommodates the industry's most generous articulation capacity and provides superior imaging for patients of diverse sizes.

Features

- 23" top slide
- 9 1/2" slide to head and 13 1/4" slide to foot
- 1,000 lb lift and 800 lb articulation weight capacity
- Removable back and leg sections for quick accessory attachment and positioning
- Power beach chair positioner for one-touch set-up
- Durable stainless steel base is easy to clean
- 243/8"-431/4" adjustable table height
- 30° lateral tilt
- 30° Trendelenburg and reverse Trendelenburg



Specifications

Length w/o head rest	69 1/2", 1,765mm
Length with head rest	82 1/8", 2,086mm
Width	19 5/8", 500mm
Width with side rails	21 7/8", 555mm
Adjustments via pendant control	
Height	24 3/8" to 43 1/4", 620 - 1,100mm
Lateral tilt (left/right)	30°
Trendelenburg	30°
Reverse Trendelenburg	30°
Power Kidney Bridge	5 3/4", 146mm Lift
Back up/down	90°/40°
Leg up/down	15/100°
Top slide/head	9 1/2", 245mm
Top slide/foot	13 1/4", 335mm
Flex/Reflex	40°/90°
Adjustments via manual control:	
Head Section	60° up/90° down
Max. permissible load	1,000 lb (450kg) lift 800 lb (360kg) Articulation